



SterilCARE™ Preventative Maintenance PartsPAX™

Consolidated offers the convenience of three PartsPAX™ for both your steam sterilizer and integral steam generator. Each PartsPAX™ includes the recommended certified components to keep your sterilizer working reliably and safely over time.



Safety PAX

Includes components that ensure your autoclave is running safely. CSS recommends this PartsPAX™ to be changed out annually regardless of number of cycles.



Performance PAX

Includes all essential components to ensure your sterilizer is running smoothly. CSS recommends this PartsPAX™ to be changed out every 300 cycles.



Complete PAX

Includes all Performance PAX and Safety PAX components plus a few other key recommended replacement parts. CSS recommends this PartsPAX™ to be used if maintenance hasn't been performed in over a year and/or the autoclave sees heavy use.



Ordering

Contact your local authorized service provider or the Consolidated parts team at parts@consteril.com.

Please provide your sterilizer serial number and an estimated number of cycles per year for your optimal SterilCARE™ PartsPAX™ or Service Plan

Each PartsPAX™ will be customized to your specific sterilizer configuration

	Safety PAX	Performance PAX	Complete PAX
Recommended Frequency	Annual Regardless of Number of Cycles	Every 300 Cycles	Annual or As Needed

Sterilizer Part/ Component

C3 Chamber Cleaner Starter Kit	✓		✓
Door Sealant	✓		✓
Internal Door Lock Assembly *	✓		✓
Door Gasket	✓		✓
Jacket Safety Valve	✓		✓
Chamber Safety Valve	✓		✓
Door Lubricant	✓		✓
Valve Assembly: Chamber, Jacket, Exhaust, Water, Gasket Air**, Gasket Exhaust**		✓	✓
Vacuum Release Assembly (Valve Repair Kit, Air Filter, Check Valve)		✓	✓
Check valves (3)			✓
Jacket Steam Trap Kit			✓
Chamber Steam Trap Kit			✓

Steam Generator Part/Component (if applicable)

Generator Safety Valve	✓		✓
Contactors	✓	✓	✓
Water Fill Valve Repair Kit		✓	✓
Sight Glass with Gasket & Valves		✓	✓
Generator Flange O-Ring		✓	✓
Low Water Cut Off Complete Assembly			✓

* For radial arm door only ** For vertical door autoclaves

SterilCARE™ Preventative Maintenance Service Plans



The Importance of Maintenance

Is a **safe** and **reliable** steam sterilizer important to your facility? In order to ensure your autoclave is running reliably and safely, Consolidated offers two simple SterilCARE™ service plans: Safety Plus and Proactive. Consolidated Sterilizer Systems recommends a safety inspection **annually** (**Safety Plus**) and a preventative maintenance inspection by an authorized factory trained service representative once every **300 cycles** (**Proactive**).

Choose the SterilCARE™ plan that best suits your facility

SterilCARE™ Plan*	Inspection Frequency	Replacement Frequency of Safety Components	Replacement Frequency of PM Components	PartsPAX Included
Safety Plus Inspection that ensures a safe sterilizer operation	Annually	Annually	As needed and determined upon inspection	Safety PAX
Proactive Planned maintenance that minimizes sterilizer downtime	Every 300 cycles	Annually	Replaced at every inspection	To be determined and based upon frequency of service

* For a customized service program please contact Consolidated to discuss your specific facility requirements.