



## SterilCARE™ Preventative Maintenance PartsPAX™

Consolidated offers the convenience of three PartsPAX™ for both your steam sterilizer and integral steam generator. Each PartsPAX™ includes the recommended certified components to keep your sterilizer working reliably and safely overtime.



### Safety PAX

Includes components that ensure your autoclave is running safely. CSS recommends this PartsPAX™ to be changed out annually regardless of number of cycles.



### Performance PAX

Includes all essential components to ensure your sterilizer is running smoothly. CSS recommends this PartsPAX™ to be changed out every 300 cycles.



### Complete PAX

Includes all Performance PAX and Safety PAX components plus a few other key recommended replacement parts. CSS recommends this PartsPAX™ to be used if maintenance hasn't been performed in over a year and/or the autoclave sees heavy use (over 1500 cycles per year).

### Each PartsPAX™ will be customized to your specific sterilizer configuration

	Safety PAX	Performance PAX	Complete PAX
Recommended Frequency	Annual Regardless of Number of Cycles	Every 300 cycles	Annual or As Needed

#### Sterilizer Part/ Component

Sterilizer Part/ Component	Safety PAX	Performance PAX	Complete PAX
C3 Chamber Cleaner Starter Kit	✓		✓
Door Sealant	✓		✓
Internal Door Lock Assembly *	✓		✓
Door Gasket	✓		✓
Jacket Safety Valve	✓		✓
Chamber Safety Valve	✓		✓
Door Lubricant	✓	✓	✓
Valve Assembly: Chamber, Jacket, Exhaust, Water, Gasket Air**, Gasket Exhaust**		✓	✓
Vacuum Release Assembly (valve repair kit, air filter, check valve)		✓	✓
Check valves (3)			✓
Jacket Steam Trap Kit			✓
Chamber Steam Trap Kit			✓

#### Steam Generator Part/Component (if applicable)

Steam Generator Part/Component (if applicable)	Safety PAX	Performance PAX	Complete PAX
Generator Safety Valve	✓		✓
Contactors	✓	✓	✓
Water Fill Valve Repair Kit		✓	✓
Sight Glass with Gasket & Valves		✓	✓
Generator Flange O-Ring			✓
Generator Circuit Board			✓
Generator Probe			✓

\* For radial arm door only \*\* For vertical door autoclaves



### Ordering

Contact your local authorized service provider or the Consolidated parts team at [parts@consteril.com](mailto:parts@consteril.com).

Please provide your sterilizer serial number and an estimated number of cycles per year for your optimal SterilCARE™ PartsPAX™ or Service Plan

## SterilCARE™ Preventative Maintenance Service Plans



### The Importance of Maintenance

Is a **safe** and **reliable** steam sterilizer important to your facility? In order to ensure your autoclave is running reliably and safely, Consolidated offers two simple SterilCARE™ service plans: Safety Plus and Proactive. Consolidated Sterilizer Systems recommends a safety inspection annually (**Safety Plus**) and a preventative maintenance inspection by an authorized factory trained service representative once every 300 cycles (**Proactive**).

### Choose the SterilCARE™ plan that best suits your facility

SterilCARE™ Plan*	Inspection Frequency	Replacement Frequency of Safety Components	Replacement Frequency of PM Components	PAX Included
<b>Safety Plus</b> Inspection that ensures a safe sterilizer operation	Annually	Annually	As needed and determined upon inspection	Safety PAX
<b>Proactive</b> Planned maintenance that minimizes sterilizer downtime	Every 300 cycles	Annually	Replaced at every inspection	To be determined and based upon frequency of service

\* For a customized service program please contact Consolidated to discuss your specific facility requirements.