



HC Complete PAX

All essential PM components for an optimal autoclaving experience. CSS recommends this Complete PAX™ to be used if maintenance hasn't been performed in over a year and/or the autoclave sees heavy use.



HC Performance PAX

Several critical components required for reliable sterilizer operation. This PartsPAX is comprised of additional OEM parts to supplement our Complete PartsPAX.



Ordering

Contact your local authorized service provider for on-site service. Or reach out to us directly for support:

617-782-6072

Parts: Ext 6 | parts@consteril.com Service: Ext 7 | service@consteril.com

Please have your sterilizer's serial # ready when you contact us for support.

SterilCARE™ Parts for Preventative Maintenance

Consolidated offers the convenience of two PartsPAX[™] for both your steam sterilizer and integral steam generator.

Each PartsPAX™ will be customized to your specific sterilizer configuration

Sterilizer Part/ Component	HC Complete PAX	HC Performance PAX
C3 Chamber Cleaner Starter Kit		
Door Gasket	✓	
Jacket Safety Valve	✓	
Chamber Safety Valve	✓	
Door Lubricant	✓	
Valve Assembly: Chamber, Jacket, Exhaust, Water, Gasket Air, Gasket Exhaust	✓	/
Vacuum Release Assembly (Valve Repair Kit, Air Filter, Check Valve)	✓	/
Check valves	✓	
Jacket Steam Trap	✓	
Chamber Steam Trap	1	
Steam Generator Part/Component	(if applicable)	
Generator Safety Valve	✓	
Contactor	✓	✓
Water Fill Valve Repair Kit	√	1

SterilCARE™ Service for Preventative Maintenance

Sight Glass with Gasket & Valves

Generator Flange 0-Ring

Low Water Cut Off



The Importance of Maintenance

To ensure your autoclave is running reliably and safely, Consolidated recommends quarterly preventative maintenance with our Proactive SterilCARE plan.

Scheduled quarterly PM visits with a factory-trained service provider will eliminate downtime and make sterilization one less thing for you to worry about.

SterilCARE™ Plan*	Inspection Frequency	Replacement Frequency of PM Components	PartsPAX Included
Proactive Planned maintenance that minimizes sterilizer downtime	2x/year	As needed and determined upon inspection	HC Complete PAX
Proactive Plus Planned maintenance, including emergency service, that minimizes sterilizer downtime.	Quarterly	Replaced at every inspection	HC Complete PAX and HC Performance PAX

^{*} For a customized service program please contact Consolidated to discuss your specific facility requirements.

